General Implants

GmbH Deutschland

Dental



EN

Experience and passion in medical manufacturing

MADE IN GERMANY



Innovative Implantology

General Implants GmbH Deutschland is a company focused on the development and production of superior implant systems that offer the utmost in technological excellence. All our systems, instruments and accessories are designed to suit the needs of contemporary implantology. Innovative strength, combined with up-to-the-minute insight into implantology advancements, guarantee the quality, safety and reliability that our partners expect.Made in Germany.

Customers are our focus

The everlasting pursue for an ideal solution to clinical and dental challenges is our daily work. Our goal is to know exactly what our customers are looking for in order to develop the perfect implant system. At General Implants, we have created a customer-orientated structure with clear communication. Customers are provided with fast, direct, informed responses. Ultimately this is just one of the reasons for our growing success.

Get to know us better
by visiting our website!


www.general-implants.com





CONTENT

Implants	6
Prosthetics	10
Drills and Drillstops	16
Tools	18
Surgical Set	20
Easy Fix	24
Fixed on Four	26



IMPLANTATE
EASY FAST **D**

SHAPE

Conical implant body adapts to complex anatomical bone structures

Tapered neck design that reduces pressure on cortical bone, ensuring better load distribution

THREAD

Variable thread for superior grip and effortless insertion


Double-thread with strong end compression for maximum stability

SURFACE

Dual Treatment Process: Combines sandblasting with biocompatible zirconium oxide and precision-controlled acid-etching

Optimized Surface Roughness: Homogeneous micro- and macro-textures designed to maximize bone cell adhesion, proliferation, and stability

Enhanced Osseointegration: Promotes faster and stronger bone-to-implant bonding, supporting long-term implant stability







EASY FAST D

Product Name	Color Code	Size	Length	Material	Thread	Connection	Product Code
Easy Fast D STERILE R	●	3.5	8.5	Titanium	Double Thread	Internal-Hex ⬡ 2.4 mm across flats	ID35085
			10.0				ID35100
			11.5				ID35115
			13.0				ID35130
			16.0				ID35160
	●	4.0	8.5				ID40085
			10.0				ID40100
			11.5				ID40115
			13.0				ID40130
			16.0				ID40160
	●	4.3	8.5				ID43085
			10.0				ID43100
			11.5				ID43115
			13.0				ID43130
			16.0				ID43160
	●	5.0	8.5				ID50085
			10.0				ID50100
			11.5				ID50115
			13.0				ID50130
			16.0				ID50160

IMPLANTATE
EASY FAST **S**

BONE SAVING

With a cylindrical body design, the Easy Fast S implant system exemplifies simplicity. It not only preserves bone substance but also promotes long-term peri-implant health

MODULARITY

The Easy Fast S Implant system has been proven to be suitable for all indications of dental implantology. The abutment system comprises all the requisite components and shares full compatibility with Easy Fast D system.

NECK



Three different implant neck designs provide an optimal biophysical surface for osseointegration at all insertion depths.





EASY FAST **S**

Product Name	Color Code	Size	Length	Connection	Material	Product Code
Easy Fast S STERILE R	●	3.0	8.5	Internal hex ⬡ 2.0 mm across flats	Titanium	IS30085
			10.0			IS30100
			11.5			IS30115
			13.0			IS30130
			16.0			IS30160
	●	3.8	8.5	Internal hex ⬡ 2.4 mm across flats		IS38085
			10.0			IS38100
			11.5			IS38115
			13.0			IS38130
			16.0			IS38160
	●	4.3	6.0			IS43060
			8.5			IS43085
			10.0			IS43100
			11.5			IS43115
			13.0			IS43130
			16.0			IS43160
	●	5.0	6.0			IS50060
			8.5			IS50085
			10.0			IS50100
			11.5			IS50115
			13.0			IS50130
			16.0			IS50160



EASY FAST **S** MTF

Product Name	Color Code	Size	Length	Connection	Material	Product Code
Easy Fast S MTF STERILE R	●	3.8	8.5	Internal hex ⬡ 2.4mm across flats	Titanium	IS38085MTF
			10.0			IS38100MTF
			11.5			IS38115MTF
			13.0			IS38130MTF
			16.0			IS38160MTF
	●	4.3	6.0			IS43060MTF
			8.5			IS43085MTF
			10.0			IS43100MTF
			11.5			IS43115MTF
			13.0			IS43130MTF
			16.0			IS43160MTF
	●	5.0	6.0			IS50060MTF
			8.5			IS50085MTF
			10.0			IS50100MTF
			11.5			IS50115MTF
			13.0			IS50130MTF
			16.0			IS50160MTF



EASY FAST **S** MF

Product Name	Color Code	Size	Length	Connection	Material	Product Code
Easy Fast S MF STERILE R	●	3.8	8.5	Internal hex ⬡ 2.4mm across flats	Titanium	IS38085MF
			10.0			IS38100MF
			11.5			IS38115MF
			13.0			IS38130MF
			16.0			IS38160MF
	●	4.3	6.0			IS43060MF
			8.5			IS43085MF
			10.0			IS43100MF
			11.5			IS43115MF
			13.0			IS43130MF
			16.0			IS43160MF
	●	5.0	6.0			IS50060MF
			8.5			IS50085MF
			10.0			IS50100MF
			11.5			IS50115MF
			13.0			IS50130MF
			16.0			IS50160MF

Notes

- Some product images are greatly magnified for better detail.
- We reserve the right to make technical changes as well as changes in shape, color and/or weight within reasonable limits.
- Not all products are available in all countries; please contact your local distributor.

PROSTHETICS

GINGIVA FORMERS






Product Name	ø	Length (L)	Material	Product Code
Gingiva former cylindrical for implant ø 3.0	3.7	2.0	Titanium	P30402-02
		3.0		P30402-03
		4.0		P30402-04
		5.0		P30402-05
		6.0		P30402-06
		7.0		P30402-07
Gingiva former cylindrical for implant ø 3.5/3.8/4.0/4.3/5.0	4.1	2.0		P44402-02
		3.0		P44402-03
		4.0		P44402-04
		5.0		P44402-05
		6.0		P44402-06
		7.0		P44402-07
Gingiva former tapered for implant ø 3.0	3.4	2.0	Titanium	P30402k-02
	3.7	3.0		P30402k-03
	4.0	4.0		P30402k-04
	4.3	5.0		P30402k-05
	4.6	6.0		P30402k-06
	4.8	7.0		P30402k-07
Gingiva former tapered for implant ø 3.5/3.8/4.0/4.3/5.0	4.8	2.0		P44402k-02
	4.9	3.0		P44402k-03
	5.0	4.0		P44402k-04
	5.2	5.0		P44402k-05
	5.3	6.0		P44402k-06
	5.5	7.0		P44402k-07

COVER SCREWS



Product Name	Thread	ø	Length	Material	Product Code
Cover screw for implant ø 3.0	M1.4	2.5	4.3	Titanium	S30460
Cover screw for implant ø 3.5/3.8/4.0/4.3/5.0	UNF1/72	3.5	5.0		S44460
Cover screw - sinus lift for implant ø 3.5/3.8/4.0/4.3/5.0	UNF1/72	6.0	5.3		S44461

PARALLEL PINS




Product Name	ø	Length	Material	Product Code
Parallel pin 0°	2.5	22	Titanium	L44650
Parallel pin 0° (4 pieces)				L44651
Parallel pin 15°				L44652
Parallel pin 15° (4 pieces)				L44653
Parallel pin 25°				L446544
Parallel pin 25° (4 pieces)				L446545


IMPRESSION POSTS











Product Name	ø	Length	Material	Product Code
Repositioning cap, closed coping	4.3	7.0	POM	L43700
Repositioning cap, closed coping, 10 pcs.	4.3	7.0		L43710
Impression post for implant ø 3.0, closed coping *	3.5	11.3	Titanium	L30512
Impression post for implant ø 3.5/3.8/4.0/4.3/5.0, closed coping *	3.5	11.6		L43512-11
Impression post for implant ø 3.5/3.8/4.0/4.3/5.0, closed coping *	4.5	11.6		L44512-11
Impression post for implant ø 3.0, open coping	3.5	14.0		L30413
Impression post for implant ø 3.5/3.8/4.0/4.3/5.0, open coping (short)	3.5	12.8		L43513-11
	4.5	12.8		L44513-11
Impression post for implant ø 3.5/3.8/4.0/4.3/5.0, open coping (long)	3.5	15.9		L43513-18
	4.5	15.9		L44513-18
Impression screw for implant ø 3.0, open coping	M1.4	23.5		S30413
Impression screw for implant ø 3.5/3.8/4.0/4.3/5.0, open coping (short)	UNF 1/72	21.4		S44413-21
Impression screw for implant ø 3.5/3.8/4.0/4.3/5.0, open coping (long)	UNF 1/72	24.8		S44413-26
Implant analog for implant ø 3.0	3.0	13.0		L30422
Implant analog for implant ø 3.5/3.8/4.0/4.3/5.0	3.8	13.0		L44422

*Screw included


PROSTHETICS

ABUTMENTS


*Screw included.



Product Name	ø	Length (L_1/L_2)	Material	Product Code
Abutment 0° * for implant ø 3.0	3.7	8.3/1.0	Titanium	P30427-01
		9.3/2.0		P30427-02
		10.3/3.0		P30427-03
		9.3/2.0	PEEK	P30427P-02
		9.3/2.0	POM	P30427Po-02
Abutment 0° (short)* for implant ø 3.5/3.8/4.0/4.3/5.0	4.5	8.6/1.0	Titanium	P44426-01
		9.6/2.0		P44426-02
		10.6/3.0		P44426-03
		9.6/2.0	PEEK	P44426P-02
		9.6/2.0	POM	P44426Po-02
Abutment 0° (long)* for implant ø 3.5/3.8/4.0/4.3/5.0	4.5	11.0/1.0	Titanium	P44427-01
		12.0/2.0		P44427-02
		13.0/3.0		P44427-03
		12.0/2.0	PEEK	P44427P-02
		12.0/2.0	POM	P44427Po-02
Abutment 15° * for implant ø 3.0	3.3	11.5/2.0	Titanium	P30432-02
	3.3	11.5/3.0		P30432-03
Abutment 15° * for implant ø 3.5/3.8/4.0/4.3/5.0	4.1	11.5/2.0		P44432-02
	4.1	11.5/3.0		P44432-03
Abutment 25° * for implant ø 3.0	3.3	11.5/2.0		P30425-02
	3.3	11.5/3.0		P30425-03
Abutment 25° * for implant ø3.5/3.8/4.0/4.3/5.0	4.1	11.5/2.0		P44425-02
	4.1	11.5/3.0		P44425-03




PREMOLAR ABUTMENTS




Product Name	ø	Length	Material	Product Code
Premolar abutment (short)* for implant ø 3.0	3.0	10.0	Titanium	P30430-10
Premolar abutment (long)* for implant ø 3.0		13.0		P30430-13
Premolar abutment (short)* for implant ø 3.5/3.8/4.0/4.3/5.0	5.5	10.0		P44430-10
Premolar abutment (long)* for implant ø 3.5/3.8/4.0/4.3/5.0		13.0		P44430-13

CAD/CAM



Product Name	ø	Length	Material	Product Code
Scan adapter* for implant ø 3.0	4.0	12.3	PEEK	P30424P-04
Scan adapter* for implant ø 3.5/3.8/4.0/4.3/5.0	4.0	12.6		P44424P-04
Scan adapter* for implant ø 3.5/3.8/4.0/4.3/5.0	6.0	12.6		P44424P-06




Product Name	ø	Length	Material	Product Code
Bonding base* for implant ø 3.0	4.0	7.6	Titanium	P30423
Bonding base* for implant ø 3.5/3.8/4.0/4.3/5.0		7.5		P44423



CASTABLE ABUTMENTS


Product Name	ø	Length	Material	Product Code
Castable abutment* for implant ø 3.0	3.7	14.0	POM	P30416
Castable abutment* for implant ø 3.5/3.8/4.0/4.3/5.0	4.6			P44416

RETENTION SCREWS



Product Name	Thread	Length	Material	Product Code
Retention screw for implant ø 3.0	M1.4	8	Titanium	S30435
Retention screw for implant ø 3.5/3.8/4.0/4.3/5.0	UNF 1/72	8.9		S44435
Retention screw (blue) for implant ø 3.0	M1.4	8		S30435b
Retention screw (blue) for implant ø 3.5/3.8/4.0/4.3/5.0	UNF 1/72	8.9		S44435b


LOCATOR ABUTMENTS



Product Name	ø	Length	Material	Product Code
Locator abutment for implant ø 3.0	3.7	2.0	Titanium	P30476-02
		3.0		P30476-03
		4.0		P30476-04
		5.0		P30476-05
		6.0		P30476-06
		7.0		P30476-07
Locator abutment for implant ø 3.5/3.8/4.0/4.3/5.0		2.0		P44476-02
		3.0		P44476-03
		4.0		P44476-04
		5.0		P44476-05
		6.0		P44476-06

PROSTHETICS

LOCATOR ANALOG



Product Name	ø	Length	Material	Product Code
Implant analog (locator)	3.7	14.0	Titanium	L44422lo

LOCATOR RETENTION ACCESSORIES



Product Name	Material	Product Code
Impression coping (locator)	Nylon	L44480



Product Name	Material	Product Code
Cap with black insert (locator)	Titanium	P44481



Product Name	Material	Product Code
Plastic insert (locator) Low Retention (4 pieces) red	Nylon	P44482-01



Product Name	Material	Product Code
Plastic insert (locator) Medium Retention (4 pieces) orange	Nylon	P44482-02



Product Name	Material	Product Code
Plastic insert (locator) High Retention (4 pieces) green	Nylon	P44482-04



Product Name	Material	Product Code
Block-Out Spacer	TPE	P44485

*Screw included

LOCATOR RETENTION ACCESSORY SET



Set Name	Content
Locator processing Set P44484	1 x Cap with black insert (locator)
	1 x Plastic insert (locator) low retention, red
	1 x Plastic insert (locator) medium retention, orange
	1 x Plastic insert (locator) low retention, green
	1 x Blockout-Out Spacer



BALL ANCHORS

Product Name	ø	Length	Material	Product Code
Transgingival ball anchor for implant ø 3.0	3.7	1.0	Titanium	P30466-01
		2.0		P30466-02
		3.0		P30466-03
		4.0		P30466-04
		5.0		P30466-05
Transgingival ball anchor for implant ø 3.5/3.8/4.0/4.3/5.0	4.5	1.0		P44466-01
		2.0		P44466-02
		3.0		P44466-03
		4.0		P44466-04
		5.0		P44466-05
		6.0		P44466-06
		7.0		P44466-07

BALL ANCHORS ACCESSORIES




Product Name	ø	Length	Material	Product Code
Impression sleeve (ball anchor)	5.95	8.0	POM	L44474
Implant analog (ball anchor)	4.5	15.8	Titanium	L44422b
Cap (ball anchor)	5.0	2.85	Titanium	P44471
Plastic insert (ball anchor) soft (4 pieces)	3.9	2.5	PEEK	P44472-01
Plastic insert (ball anchor) reinforced retention (4 pieces)	3.9	2.5	PEEK	P44472-02



BALL ANCHOR SET


Set Name	Content
Ball anchor set P44473	Cap (ball anchor)
	Plastic insert (ball anchor) soft
	Plastic insert (ball anchor) reinforced retention

DRILLS AND DRILLSTOPS



SURGICAL DRILLS (PREPARATION)

Product Name	Ø	Length	Material	Product Code
Round drill	2.3	8	Stainless steel	B2140
Triangular drill	1.8	13		B2240
Pilot drill	2.0	16		B2040




SURGICAL DRILLS (EASY FAST S)

Product Name	Color Code	Ø	Length	Material	Product Code
Regular drill (EFS ø 3.0)	●	2.5	16	Stainless steel	B2530
Regular drill (EFS ø 3.8)	●	3.2			B3238
Regular drill (EFS ø 4.3)	●	3.7			B3743
Regular drill (EFS ø 5.0)	●	4.4			B4450



TAPS (EASY FAST S)


Product Name	Color Code	Ø	Length	Material	Product Code
Tap (EFS ø 3.0)	●	3.0	4.6	Stainless steel	T90300
Tap (EFS ø 3.8)	●	3.8			T90380
Tap (EFS ø 4.3)	●	4.3			T90430
Tap (EFS ø 5.0)	●	5.0			T90500





SURGICAL DRILLS (EASY FAST D)

Product Name	Color Code	Ø	Length	Material	Product Code
Regular drill (EFD ø 3.5)	●	3.2	16	Stainless steel	B353227
Regular drill (EFD ø 4.0)	●	3.7	16		B403730
Regular drill (EFD ø 4.3)	●	3.9	16		B433933
Regular drill (EFD ø 5.0)	●	4.6	16		B504640
Cortical drill (EFD ø 3.5)	●	3.6	5		B353629
Cortical drill (EFD ø 4.0)	●	3.9	5		B403932
Cortical drill (EFD ø 4.3)	●	4.2	5		B434235
Cortical drill (EFD ø 5.0)	●	4.9	5		B504942



DRILL STOPS (EASY FAST S)


Product Name	Color Code	Ø	Length(L_1/L_2)	Material	Product Code
Drill stop for pilot drill	-	2.0	3.0/13.0	Titanium	T2040-030
			4.5/11.5		T2040-045
			6.0/10.0		T2040-060
			7.5/8.5		T2040-075
			9.7/6.3		T2040-118
Drill stop for implant ø 3.0	●	2.5	3.0/13.0		T2530-030
			4.5/11.5		T2530-045
			6.0/10.0		T2530-060
			7.5/8.5		T2530-075
Drill stop for implant ø 3.8	●	3.2	3.0/13.0		T3238-030
			4.5/11.5		T3238-045
			6.0/10.0		T3238-060
			7.5/8.5		T3238-075
Drill stop for implant ø 4.3	●	3.7	3.0/13.0		T3743-030
			4.5/11.5		T3743-045
			6.0/10.0		T3743-060
			7.5/8.5		T3743-075
			9.7/6.3		T3743-118
Drill stop for implant ø 5.0	●	4.4	3.0/13.0		T4450-030
			4.5/11.5		T4450-045
			6.0/10.0		T4450-060
			7.5/8.5		T4450-075
			9.7/6.3		T4450-118

DRILL STOPS (EASY FAST D)

Product Name	Color Code	Ø	Length(L_1/L_2)	Material	Product Code
Drill stop for implant ø 3.5	●	3.2	3.0/13.0	Titanium	T353227-030
			4.5/11.5		T353227-045
			6.0/10.0		T353227-060
			7.5/8.5		T353227-075
Drill stop for implant ø 4.0	●	3.7	3.0/13.0		T403730-030
			4.5/11.5		T403730-045
			6.0/10.0		T403730-060
			7.5/8.5		T403730-075
Drill stop for implant ø 4.3	●	3.9	3.0/13.0		T433933-030
			4.5/11.5		T433933-045
			6.0/10.0		T433933-060
			7.5/8.5		T433933-075
Drill stop for implant ø 5.0	●	4.6	3.0/13.0		T504640-030
			4.5/11.5		T504640-045
			6.0/10.0		T504640-060
			7.5/8.5		T504640-075


TOOLS

SCREWDRIVERS (IMPLANTS)




Product Name	⬡	Length(L_{total}/L)	Material	Product Code
Screwdriver for implant ø 3.0 (manual)	2.0	16/7	Stainless steel	T30563-07
		26/17		T30563-17
		36/27		T30563-27
Screwdriver for implant ø 3.5/3.8/4.0/4.3/5.0 (manual)	2.4	16/7		T44563-07
		26/17		T44563-17
		36/27		T44563-27
Screwdriver for implant ø 3.0 (contra-angle)	2.0	26/10	Stainless steel	T30564-10
		36/20		T30564-20
Screwdriver for implant ø 3.5/3.8/4.0/4.3/5.0 (contra-angle)	2.4	26/10		T44564-10
		36/20		T44564-20

SCREWDRIVERS (LOCATOR ABUTMENT)



Product Name	⬡	Length(L_{total}/L)	Material	Product Code
Screwdriver for locator abutment (manual)	2.0	16/7	Stainless steel	T20563

SCREWDRIVER (RETENTION SCREW)




Product Name	⬡	Length(L_{total}/L)	Material	Product Code
Screwdriver for retention screw (manual)	1.2	19.0/10.0	Stainless steel	T40565-10
		24.0/15.0		T40565-15
		34.0/25.0		T40565-25
Self-Handled Screwdriver for retention screw (manual)	1.2	21,1/10.0	Stainless steel	T40566-10
		26.1/15.0		T40566-15
		31.1/20.0		T40566-20
Screwdriver for retention screw (contra-angle)	1.2	23.5/10.0	Stainless steel	T40567-10
		33.5/20.0		T40567-20

MUCOSAL PUNCH



Product Name	ø	Length	Material	Product Code
Mucosal punch	3.0	21.0	Stainless steel	T30533-03
	4.0			T30533-04
	5.0			T30533-05



MISCELLANEOUS TOOLS

Product Name	Material	Product Code
Insertion tool (pin head and tube head) (ball anchor)	Stainless steel	T44470
Measuring probe	Stainless steel	T44535
Drill shaft extension	Stainless steel	T44536
Contra-angle to quad Adapter for manual tools	Stainless steel	T44579
Hand wheel	Stainless steel	T4018
Handle	Stainless steel	T44130
Ratchet	Stainless steel	T44700
Torque ratchet (10-40 Ncm)	Stainless steel	T44800
Locator tool	Stainless steel	T441496

Notes

- Some product images are greatly magnified for better detail.
- We reserve the right to make technical changes as well as changes in shape, color and/or weight within reasonable limits.
- Not all products are available in all countries, please contact your local distributor.

One Surgical-Set for all Implants

Product Name	Product Code
Surgical set PROFESSIONAL EASY FAST S/D	EF2000



EASY FAST S/D INSTRUMENTS SET

	Product Name	Product Code
1	Pilot drill	B2040
2	Triangular drill	B2240
3	Round drill	B2140
4	Regular drill (EFS ø 3.0)	B2530
5	Tap (EFS ø 3.0)	T90300
6	Regular drill (EFS ø 3.8)	B3238
7	Tap (EFS ø 3.8)	T90380
8	Regular drill (EFS ø 4.3)	B3743
9	Tap (EFS ø 4.3)	T90430
10	Regular drill (EFS ø 5.0)	B4450
11	Tap (EFS ø 5.0)	T90500
12	Regular drill (EFD ø 3.5)	B353227
13	Cortical drill (EFD ø 3.5)	B353629
14	Regular drill (EFD ø 4.3)	B433933
15	Cortical drill (EFD ø 4.3)	B434235
16	Regular drill (EFD ø 5.0)	B504640
17	Cortical drill (EFD ø 5.0)	B504942
18	Screwdriver for Implant ø 3.0 (manual)	T30563-17
19	Screwdriver for Implant ø 3.0 (manual)	T30564-10
20	Screwdriver for Implant ø 3.0 (manual)	T30564-20
21	Screwdriver for Implant ø 3.5/3.8/4.3/5.0 (manual)	T44563-07
22	Screwdriver for Implant ø 3.5/3.8/4.3/5.0 (manual)	T44563-17
23	Screwdriver for Implant ø 3.5/3.8/4.3/5.0 (contra angle)	T44564-10
24	Screwdriver for Implant ø 3.5/3.8/4.3/5.0 (contra angle)	T44564-20
25	Screwdriver for Locator abutment (manual)	T20563
26	Screwdriver for retention screw (manual)	T40565-10
27	Screwdriver for retention screw (manual)	T40565-25
28	Screwdriver for retention screw (manual)	T40566-15
29	Screwdriver for retention screw (contra angle)	T40567-10
30	Handle	T44130
31	Hand wheel	T4018
32	Torque ratchet (10-40 Ncm)	T44800
33	EASY FAST S/D Tray, empty	EF2000L

*Tool for EFD 4mm is not included

EASY FAST **D**

Compact and Modern Surgical-Set for EASY FAST **D**

Product Name	Product Code
Surgical set PROFESSIONAL EASY FAST D	EF2000-02



EASY FAST **D** INSTRUMENTS SET

	Product Name	Product Code
1	Pilot drill	B2040
2	Triangular drill	B2240
3	Round drill	B2140
4	Regular drill (EFD ø 3.5)	B353227
5	Cortical drill (EFD ø 3.5)	B353629
6	Regular drill (EFD ø 4.0)	B403730
7	Cortical drill (EFD ø 4.0)	B403932
8	Regular drill (EFD ø 4.3)	B433933
9	Cortical drill (EFD ø 4.3)	B434235
10	Regular drill (EFD ø 5.0)	B504640
11	Cortical drill (EFD ø 5.0)	B504942
12	Screwdriver for Implant ø 3.5/3.8/4.0/4.3/5.0 (manual)	T44563-07
13	Screwdriver for Implant ø 3.5/3.8/4.0/4.3/5.0 (manual)	T44563-17
14	Screwdriver for retention screw (manual)	T40565-10
15	Screwdriver for retention screw (manual)	T40565-25
16	Handle	T44130
17	Hand wheel	T4018
18	Screwdriver for retention screw (manual)	T40566-15
19	Screwdriver for Implant ø 3.5/3.8/4.0/4.3/5.0 (contra angle)	T44564-10
20	Screwdriver for Implant ø 3.5/3.8/4.0/4.3/5.0 (contra angle)	T44564-20
21	Screwdriver for retention screw (contra angle)	T40567-10
22	Torque manual (10-40 Ncm)	T44800
23	Screwdriver for Locator abutment (manual)	T20563
24	Regular drill Ø 3,0	B2530
25	EASY FAST D Tray, empty	EFD2000-02L


EASY FIX



Easy Fix facilitates immediate prosthetic treatment!

EASY FIX

The ready-made, Easy Fix bar connector system permits immediate prosthetic treatment. The Easy Fix concept can be adapted to any implant system. This system completely forgoes the use of moulded and milled connector bars. The bars can be screwed together using spacer sleeves available in different sizes. The riders are made of polycarbonate, and can be polymerised into the prosthesis using a lining.



EASY FIX BAR CONNECTOR COMPONENTS


Product Name	Ø	Length	Material	Product Code
Tube	2.0	20.9	Titanium	EFX1010
Pin	2.5	19.4	Titanium	EFX1020
Double bar	2.0/2.5	35.3	Titanium	EFX1030
Riders (3 pieces)	6.0	16.0	Polycarbo-nate	EFX1040



SPACER SLEEVES

Product Name	Ø	Length	Material	Product Code
Spacer sleeve Easy Fast S* for implant ø 3.0 (4 pieces)	4.2	1.5	Titanium	EFX3001-015
		2.5		EFX3001-025
		3.5		EFX3001-035
Spacer sleeve Easy Fast S/D* for implant ø 3.5/3.8/4.0/4.3/5.0 (4 pieces)		1.5		EFX4401-015
		2.5		EFX4401-025
		3.5		EFX4401-035

*Screw included



Product Name	Material	Product Code
Bending Tools	Stainless Steel	EFX1050

BAR CONNECTOR SETS

Size	1.5	2.5	3.5	1.5	2.5	3.5
For	*Easy Fast S ø 3.0*	*Easy Fast S ø 3.0*	*Easy Fast S ø 3.0*	*Easy Fast S/D ø 3.5-5.0*	*Easy Fast S/D ø 3.5-5.0*	*Easy Fast S/D ø 3.5-5.0*
Set No.	**EFX1300-015**	**EFX1300-025**	**EFX1300-035**	**EFX1440-015**	**EFX1440-025**	**EFX1440-035**
Tube	EFX1010	EFX1010	EFX1010	EFX1010	EFX1010	EFX1010
Pin	EFX1020	EFX1020	EFX1020	EFX1020	EFX1020	EFX1020
Double Bar x2	EFX1030	EFX1030	EFX1030	EFX1030	EFX1030	EFX1030
Riders (3 pieces)	EFX1040	EFX1040	EFX1040	EFX1040	EFX1040	EFX1040
Spacer Sleeve (4 pieces)	EFX3001-015	EFX3001-025	EFX3001-035	EFX4401-015	EFX4401-025	EFX4401-035
Bending tools	EFX1050	EFX1050	EFX1050	EFX1050	EFX1050	EFX1050

Notes

- Some product images are greatly magnified for better detail.
- We reserve the right to make technical changes as well as changes in shape, color and/or weight within reasonable limits.
- Not all products are available in all countries; please contact your local distributor.

Fixed on Four

FIXED PROSTHETICS WITH GENERAL IMPLANTS


Our Fixed on Four dental implant solutions offer an efficient approach to full-arch restoration, providing patients with a stable, and aesthetically pleasing alternative to traditional dentures.


Utilizing just four strategically placed implants, this method maximizes support and minimizes the need for extensive bone grafting, making it an ideal choice for patients with limited bone density.


The Fixed on Four system ensures faster recovery times, immediate function, and a natural-looking smile, all while delivering exceptional long-term results.


Explore our range of products designed to deliver comfort, durability, and confidence to your patients
















FIXED ON FOUR

PRODUCT NAME	Ø	LENGTH	MATERIAL	PRODUCT CODE
Abutment 0°	5.6	1.0	Titanium	FOF1001-01
		2.0		FOF1001-02
		3.0		FOF1001-03
Abutment 18,5° *	5.6	1.0	Titanium	FOF1017-01
		2.0		FOF1017-02
		3.0		FOF1017-03
Abutment 37° *	5.6	1.0	Titanium	FOF1035-01
		2.0		FOF1035-02
		3.0		FOF1035-03
Prosthetic cap *	5.6	18.0	Titanium	FOF1002
			POM	FOF1002Po
Impression cap for open impression	5.5	8.35	Titanium	FOF1003
Implant analog	5.6	15.2	Titanium	FOF1004
Screwdriver for Abutment 0° and Implant analog	SW 2.4	27.0	Stainless steel	T34205-18

*Screws included

General Implants GmbH Deutschland

Eisenbahnstrasse 100
78573 Wurmlingen
Germany

Telephone +49 (0) 7461 9336 - 0
info@general-implants.com
www.general-implants.com

CE 0483

General Implants
GmbH Deutschland

General Implants GmbH Deutschland / All rights reserved/ GI Dental Catalogue_EN_04-2025_Rev2.1